

Colorimetric Plasminogen Activation Assay

Stack Lab Protocol

(see next page for calculations)

Abbreviations:

VLKpNA - Val-Leu-Lys-para-nitroanilide (Sigma V-0882)

NA - nitrophenolate anion

Pg - plasminogen

Pm - plasmin

PA - plasminogen activator

K_m - Michaelis constant

k_{cat} - catalytic rate constant

To assay PA activity in conditioned medium using a 96 well plate microtiter plate reader:

1. Turn on plate reader before starting. On our reader, choose KINETICS405. This monitors the change in absorbance at 405 nm over a 2 hr time course at 37 C. [Turn on incubator]. Assay each sample in triplicate.

2. Reaction:

120 ul of 20 mM Hepes pH 7.4

17.5 ul of 3 uM Pg

18 ul of 3 mM VLK-pNA

Start reaction by addition of 20 ul conditioned medium sample.

Monitor A405/time.

3. To calculate mol Pm/min, see below.

To calculate the velocity of plasminogen activation
(derived by Wohl et al., (1980) J. Biol. Chem 255:2005-2013)


The reaction being measured is:




This is analyzed using a secondary reaction:



When Pg is activated by PA, cleavage of VLK-pNA is monitored by absorbance at 405 nm. If you plot absorbance vs time, a curved line results:



This is because Pm is generated continuously, so the amount of Pm in the secondary reaction is not constant. To obtain linear data, plot A_{405} vs time^2 :



This gives a straight line with a slope of $b/2$. This slope can now be used to calculate the initial velocity of Pm formation according to the equation:

$$v = \frac{b (1 + K_m/S_o)}{\epsilon k_{cat}}$$

v = initial velocity of Pg activation

K_m = Michaelis constant for hydrolysis of VLK-pNA by Pm
 $= 3 \times 10^{-4} M$

k_{cat} = catalytic rate constant for hydrolysis of VLK-pNA by plasmin.
 This constant is empirically determined in 20 mM Hepes, pH 7.4 to be $1.35 \times 10^3 M/min \cdot mol Pm$. [Note - this number may vary if your buffer is very different]

ϵ = extinction coefficient of nitrophenolate anion which absorbs at 405 nm
 $= 1 \times 10^4 M^{-1}$ [Erlanger et al., (1961) Arch. Biochem. Biophys. 95:271-278]

S_o = initial concentration of VLKpNA (0.3 mM in our assays)
 $(3 \times 10^{-4} M)$

So for example, if we plot Abs vs time², the slope of the line in A/min^2
 $= b/2 =$ (for example) 1.01×10^{-4} ; then $b = 2.02 \times 10^{-4}$

Plugging this into the velocity equation:

$$v = \frac{(2.02 \times 10^{-4} \text{ abs/min}^2) [1 + (3 \times 10^{-4} M / 3 \times 10^{-4} M)]}{(1 \times 10^4 M^{-1}) (1.32 \times 10^3 M \text{ min}^{-1} \text{ mol Pm}^{-1})}$$

$$= 3 \times 10^{-11} \text{ mol Pm/min}$$

Converting microtitre plate reader data

Our Molecular Devices microtitre plate reader gives data in milliO.D./min

1. $x \text{ mO.D./min} \times 1 \text{ O.D./1000 mO.D.} = x' \text{ O.D./min}$
2. Run time - example 120 min
 $x' \text{ O.D./min} \times 120 \text{ min} = x'' \text{ O.D.}$
3. convert to O.D./min^2 :

$$x'' \text{ O.D./120}^2\text{min} = x''' \text{ O.D./min}^2$$

Plug this into equation above to get mol Pm/min. [If using our standard assay conditions of 0.3 mM VLK-pNA and 3 mM Pg in 20 mM Hepes, pH 7.4, can just multiply by 1.5×10^{-7})

Example: Plate reader gives you 10 mO.D./min for a 2 hr run

1. $10 \text{ mO.D./min} \times 1 \text{ O.D./1000 mO.D.} = 0.01 \text{ O.D./min}$
2. $0.01 \text{ O.D./min} \times 120 \text{ min} = 1.2 \text{ O.D.}$
3. $1.2 \text{ O.D./120}^2\text{min}^2 = 8.3 \times 10^{-5} \text{ O.D./min}^2$

Plug into equation and result is $1.25 \times 10^{-11} \text{ mol Pm/min}$